

NGS Analysis and Galaxy

Part 2


Kathryn Crouch
kathryn.crouch@glasgow.ac.uk

Outline

- Background
 - How does sequencing work?
 - What does the data look like?
- Pre-processing for Analysis
 - QC and trimming
 - Aligning to a reference genome
- Whole-genome sequencing
 - Calling SNPs
 - CNVs

How Does Sequencing Work?

- Carry out replication under controlled conditions
- Artificially slow down the reaction to see the order in which bases are incorporated




How Does Sequencing Work?

Sanger Sequencing (1975)

Uses modified nucleotides (ddNTPs) that cannot be extended

Each ddNTP is labelled with a different dye so you can see the order in which they are incorporated

Long reads and low error rate, but low throughput




How Does Sequencing Work?

Illumina Sequencing (2005)


DNA fragments are
adaptor-ligated and attached to
a flow cell

PCR is carried out in situ to form
clusters


Sequencing can be carried out
on millions of clusters
simultaneously




illumina



How Does Sequencing Work?



How Does Sequencing Work?




How Does Sequencing Work?

Illumina Sequencing

Blocked and labelled nucleotides are added


1 nucleotide is incorporated and an image is taken of the array

Label and block are removed and cycle repeats




How Does Sequencing Work?

T G C T A C G A T



C A A T A G A C G


How Does Sequencing Work?



Paired End Illumina Sequencing

A short read is sequenced from
each end of each fragment

Blue and black are adaptors




Paired End Sequencing

A short read is sequenced from each end of each fragment

Blue and black are adaptors

Insert or fragment size is the distance between adaptors




Paired End Sequencing

A short read is sequenced from each end of each fragment

Blue and black are adaptors

Insert or fragment size is the distance between adaptors

The region between the reads is not sequenced - it is covered by other fragments



Paired End Sequencing

A short read is sequenced from each end of each fragment


Blue and black are adaptors

Insert or fragment size is the distance between adaptors

The region between the reads is not sequenced - it is covered by other fragments

Reads that map in the correct orientation and the expected distance apart are “concordant” or “proper pairs”

Concordant alignments are prioritised



File Formats: FASTQ

```
@NS500205:127:HW272BGXX:1:11101:7788:1040 1:N:0:TCTCGCGC+AGGCTATA
```

For paired-end reads you will have two files

4 lines per read

- Line 1 is a unique header (this will be shared between the pairs)

File Formats: FASTQ

```
@NS500205:127:HW272BGXX:1:11101:7788:1040 1:N:0:TCTCGCGC+AGGCTATA  
CCTTGNTCCGTCATATTTTTTAGCATTGCAATGACGCTAAGTCCCGATTGACGCGCACGTGCTCACCCGGTTTCC
```

For paired-end reads you will have two files

4 lines per read

- Line 1 is a unique header (this will be shared between the pairs)
- Line 2 is the sequence of the read

File Formats: FASTQ

```
@NS500205:127:HW272BGXX:1:11101:7788:1040 1:N:0:TCTCGCGC+AGGCTATA
CCTTGNTCCGTCATATTTTTTAGCATTGCAATGACGCTAAGTCCCGATTGACGCGCACGTGCTCACCCGGTTTCC
+
AAAAA#EEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEEE
```

For paired-end reads you will have two files

4 lines per read

- Line 1 is a unique header (this will be shared between the pairs)
- Line 2 is the sequence of the read
- Line 4 is the quality for each base
 - Quality is encoded using ASCII
 - <http://www.asciitable.com/>

Phred quality scores are logarithmically linked to error probabilities

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%
60	1 in 1,000,000	99.9999%

Alignment-based Analysis

Quality
Control



Read
Trimming



Alignment

Analysis of NGS Sequencing

Quality
Control




Read
Trimming




Alignment

- FASTQC
 - <https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>
 - Overall sequencing quality
 - GC content
 - N content
 - Read length distribution
 - Over-represented sequences
 - Adaptor content
- Output is an html file that can be opened in a web browser

✓ Per base sequence quality




! Per sequence GC content




! Overrepresented sequences

Sequence	Count	Percentage	Possible Source
ACAAGTGTGTAACATTAATTGCAAGTTTGCAACGCTGTTCTTTAGTGTT	70896	0.12562741276052788	No Hit


Analysis of NGS Sequencing




- Trimmomatic
<https://github.com/usadellab/Trimmomatic>
- Sequence quality tends to decrease towards the 3' end of reads - these can affect mapping
- Sickle will:
 - Remove poor quality reads from the 3' end of each read
 - Check for reads that are too short and discard them
 - Check that all reads still have a pair and discard those that don't
- Some trimming tools can also remove adaptors




Analysis of NGS Sequencing




- Trimmomatic
<https://github.com/usadellab/Trimmomatic>
- Sequence quality tends to decrease towards the 3' end of reads - these can affect mapping
- Sickle will:
 - Remove poor quality reads from the 3' end of each read
 - Check for reads that are too short and discard them
 - Check that all reads still have a pair and discard those that don't
- Some trimming tools can also remove adaptors



Analysis of NGS Sequencing



- Trimmomatic
<https://github.com/usadellab/Trimmomatic>
- Sequence quality tends to decrease towards the 3' end of reads - these can affect mapping
- Sickle will:
 - Remove poor quality reads from the 3' end of each read
 - Check for reads that are too short and discard them
 - Check that all reads still have a pair and discard those that don't
- Some trimming tools can also remove adaptors



Analysis of NGS Sequencing

Quality
Control



Read
Trimming



Alignment

What is an alignment?

Two sequences:

ATTGAAAGCTA
GAAATGAAAAGG

How would you align one to the other?

Analysis of NGS Sequencing

Quality
Control



Read
Trimming



Alignment

What is an alignment?

Two sequences:

ATTGAAAGCTA
GAAATGAAAAGG

How would you align one to the other?

```
--ATTGAAA-GCTA
  | | | | | |
GAAATGAAAAGG
```

```
ATTGAAA-GCTA---
  | | | | |
---GAAATGAAAAGG
```

Which one is better??

Analysis of NGS Sequencing

Quality
Control



Read
Trimming



Alignment

What is an alignment?

Two sequences:

ATTGAAAGCTA
GAAATGAAAAGG

How would you align one to the other?

```
--ATTGAAA-GCTA
  | | | | | |
GAAATGAAAAGG--
```

```
ATTGAAA-GCTA---
  | | | | |
---GAAATGAAAAGG
```

Which one is better??

Alignment scoring:

- 1 for a match
- -1 for a mismatch
- -2 for a gap

Now alignment 1 scores -4 and alignment 2 scores -10 so we would choose alignment 1

Alignment algorithms are designed to align data in a reasonable time on a standard computer


Heuristic - not exhaustive, but “good enough”

Improving Alignments with Paired End Reads


Paired end reads can resolve alignments in repetitive regions

A read that could map in multiple locations due to repetitions in the genome can be located accurately by inference from the position of its pair

In this case, we know that the insert size is ~400 bp, so we can infer that the first alignment is more likely to be correct.



Visualising Alignments (Artemis)



Alignment Tools

Quality
Control



Read
Trimming



Alignment

Mapping Tools for DNA-seq data

- BWA
 - <https://github.com/lh3/bwa>
- Bowtie2
 - <https://bowtie-bio.sourceforge.net/bowtie2/index.shtml>

Whole Genome Sequencing (DNA-seq, WGS)

What Can We Discover From Aligned Reads?

- Where and how is our sample different from the reference?
 - Discovery of SNVs and Indels
- Coverage
 - Discovery of copy number variations
 - Discovery of regions of high variability

Finding SNVs

Quality
Control




Read
Trimming



Alignment



SNP Calling



Finding Variants

If we load alignments into a genome viewer, we can see variants

How do we find them globally? How do we assess them?

Analysis of NGS Sequencing

Quality
Control




Read
Trimming



Alignment



SNP Calling



Blue lines are reads aligned against a reference (black). Red dots indicate individual bases where a base in a read differs from the reference.

A: Most reads differ from the reference -> homozygous SNP

B and C: Roughly 50% of reads differ from the reference -> potential heterozygous SNP

D: Only one base differs from the reference -> probably a sequencing error

E: All bases the same as the reference

Analysis of NGS Sequencing

Quality
Control




Read
Trimming



Alignment



SNP Calling



Things to think about:

- What happens if your sample is not a clone?
- What happens if your sequencing depth is low?

Analysis of NGS Sequencing

Quality
Control




Read
Trimming



Alignment



SNP Calling



Freebayes <https://github.com/freebayes/freebayes>

Automated tool to call SNPs

You may also come across other tools including GATK and BCFTools.

What Else Can We Find Out?

Quality
Control



Read
Trimming



Alignment




SNP Calling


Coverage

- Expect coverage to be even across the genome
- In reality, we see variation associated with:
 - GC content
 - Repetitive or highly variable regions
 - Large scale insertions and deletions
- Note that doing alignments and examining coverage is the basis of RNAseq and ChIPseq analysis too!


Coverage and Copy Number Variations



Copy Number Variations



Copy Number Variations




Reads sampled from both copies of the duplication will all map to the single copy in the reference, increasing the coverage at this locus




Reference genome without duplication

Variable
coverage in
repetitive
regions

43.8




Segmental Duplication in *C. albicans*



Segmental tetraploidy on the right arm of chromosome 2 in a clinical sample

Galaxy



Accessing Data

- Workflows failed because we ran out of disc space - you didn't do anything wrong!
- We will explore the output of workflows that Eve pre-ran
- We will do this in the live Galaxy site NOT the workshop Galaxy site!